1、SEND AI LANGUAGE INSTRUCTIONS

1,2,3: 语音状态 0:无意义,1:开始(询问),2:结束,3:可传输 / Voice status 0: Meaningless, 1: Start (Question), 2: End, 3: Transmissible

SEND AI VOICE-TO-TEXT TEXT

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SEND AI ANSWER TEXT

text: app 回复的文本 / The text replied by the app

1,2,3
SEND AI LANGUAGE INSTRUCTIONS
text
SEND AI VOICE-TO-TEXT TEXT
text
SEND AI ANSWER TEXT

2、SEND AI TO DISPLAY HEART RATE,BLOOD OXYGEN, ACTIVITY AND OTHER DATA UI: 视图模板
/View template

title: 标题 /for example: Average Heart Rate

value: 值 / for example: 100

unit : 单位 /for example: bpm

footer: 下标 /for example: Yesterday

action_time : 时间戳 (年月日时分秒) /Timestamp (year, month, day, hour, minute, second)

title
value
unit
footer
2000-01-01 00:00:00
SEND AI TO DISPLAY HEART RATE, BLOOD OXYGEN, ACTIVITY AND OTHER DATA UI

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4、SEND AI ACTION BEHAVIOR: 动作模板/Action Template

- directive 0: 设置低心率警报; Set_Low_HR_Alert
 1: 设置高心率警报; Set_High_HR_Alert
 2: 设置低血氧饱和度警报; Set_Low_SpO2_Alert
 3: 设置高压力警报; High_Stress_Alert
 4: 设定睡眠目标; Set_Sleep_Goal
 5: 设置步数目标; Set_Steps_Goal
 6: 设定距离目标; Set_Distance_Goal
 7: 设定活动小时目标; Set_Active_Hours_Goal
 8: 设定卡路里目标; Set_Calories_Goal
 9: 设定站立目标; Set_Standing_Goal
 10: 上一经期开始日期; Log_Last_Period_Date
 11: 经期长度; Log_Period_Length
 12: 周期长度; Log_Cycle_Length

threshold_value: 详细的值, 包括目标值/Detailed values, including target values

time_date: 日期(年月日)暂时只有记录上一经期开始日期需要/The date (year, month, day) is currently only required when the start date of the record is specified.

unit: 单位/unit

action_time: 时间戳（年月日时分秒）/Timestamp (year, month, day, hour, minute, second)
参考表格/Reference table:

Scenario 场景	type 类型	target_metric 目标指标	condition 条件	threshold_value 阈值	unit 单位	页面
Set Low HR Alert 设置低 HR 警报	"alert" "警报"	"heart_rate" "心率"	"Lower Limit" "下限值"	"50"	"BPM" "BPM"	High HR alert
Set High HR Alert 设置高HR警报	"alert" "警报"	"heart_rate" "心率"	"Upper Limit" "上限"	"120"	"BPM" "BPM"	Low HR alert
Set Low SpO2 Alert 设置低 SpO2 警报	"alert" "警报"	"spo2" "血氧饱和度"	"Lower Limit" "下限值"	"90"	"%" "%"	Low SpO2 alerts
Set High Stress Alert 设置高压力警报	"alert" "警报"	"stress" "压力"	"Upper Limit" "上限"	"80"	"%" "%"	Relax reminder
Set Sleep Goal 设置睡眠目标	"goal" "目标"	"sleep" "睡眠"	"Target" "目标"	"8"	"hours" "小时"	睡眠
Set Steps Goal 设置步数目标	"goal" "目标"	"steps" "步数"	"Target" "目标"	"10000"	"steps" "步数"	活动
Set Distance Goal 设定距离目标	"goal" "目标"	"distance" "距离"	"Target" "目标"	"5"	"km" "公里"	
Set Calories Goal 设置卡路里目标	"goal" "目标"	"calories" "卡路里"	"Target" "目标"	"500"	"kcal" "千卡"	
Set Standing Goal 设定站立目标	"goal" "目标"	"standing_hours" "站立小时数"	"Target" "目标"	"12"	"hours" "小时"	
Log Last Period Date 上次记录日期	"log" "日志"	"menstrual_cycle" "月经周期"	"Last Period Start" "上一周期开始时间"	"2025-06-10"	"date" "日期"	经期对应时间设置界面
Log Period Length 日志周期长度	"log" "日志"	"menstrual_cycle" "月经周期"	"Period Length" "周期长度"	"5"	"days" "天"	
Log Cycle Length 日志周期长度	"log" "日志"	"menstrual_cycle" "月经周期"	"Cycle Length" "周期长度"	"28"	"days" "天"	

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contacts_number: 输入联系人号码/Enter the contact person's phone number

workout_type: 输入运动类型/Input movement type

Duration: 输入计时器秒数, 单位秒/Enter the timer duration in seconds, with the unit being seconds.

music_command: 音乐控制 (0:播放,1:暂停,2:下一首,3:上一首)/Music Control (0: Play, 1: Pause, 2: Next, 3: Prev)

Activate the
function template

scenario
2000-01-01 00:00:00
contactsNumber
workoutType
duration
musicCommand
SEND AI TO OPEN DEVICE RELATED SETTINGS